


Restalyst

Abundant Life Through Restfulness













AGENDA

- Introduction
 - Company Profile
 - Milestones
 - Products & Quality Assurance
- Products REAAD™ Series
 - − NPC-REAAD™
 - − GC-REAAD™
 - − HCC-REAAD™
- Products COVID-19 Tests
 - RAPIT COVID-19 IgM/IgG Test Kit
 - − COVID19N-REAAD™
- Services & Collaborations



GC-REAAD™ ITIH3 ELISA



INTRODUCTION

COMPANY PROFILE

Restalyst

- a member of the Reste Group
- privately-owned biomedical clinical diagnostic company from Singapore with a global presence.



We are committed to:

- Provide innovative, reliable and clinically-proven diagnostic solutions
- Our mission of harnessing technological advancements for the improvement of disease diagnosis and management of patient health

COMPANY PROFILE

- Fully integrated company established in 2007
- Capable of generating the full product pipeline including;
 - Scientific research
 - In-house developmental work of patented technology
 - Product manufacturing
- Research & Manufacturing facilities are EN ISO13485:2016 certified



MILESTONES

- 2007: Restalyst Pte Ltd was founded
- 2009: R&D and production facilities achieved ISO13485 certification
- 2011: NPC-REAAD received CE Mark
- 2012: NPC-REAAD listed on Singapore Medical Device Register
- 2015: Launched GC-REAAD & received CE mark
- 2016: Launched HCC-REAAD & received CE mark
- 2017: GC-REAAD listed on Singapore Medical Device Register
- 2019: HCC-REAAD listed on Singapore Medical Device Register
- 2020: Launched RAPIT COVID-19 IgM/IgG
 Test Kit & received CE & HSA PA
- 2020: Launched COVID19N-REAAD & received CE & HSA PA



PRODUCTS & QUALITY ASSURANCE

- Design & manufacturing of IVD medical devices for testing of markers for cancer & infectious diseases
 - − NPC-REAAD™
 - GC-REAAD™
 - − HCC-REAAD™
 - RAPIT COVID-19 IgM/IgG Test Kit
- ISO13485:2016, EN ISO13485:2016 (since 2009)
- Singapore Medical Device Register Listing
- CE mark
 - Certificate No: CE/SGP/2014/12/01





PRODUCTS - REAADTM SERIES

REAADTM Series

- A range of innovative diagnostic solutions for early detection of oncology diseases
- Employ the use patented or proprietary protein biomarkers
- Clinical-based approach to generate evidence-based results
- Intended for the early diagnosis for diseases



- Products span from specialised manual testing to full integration onto automated routine clinical laboratory testing solutions
- Our clinically-proven diagnostic solutions includes:
 - NPC-REAAD™ EBV EA IgA ELISA
 - GC-REAAD™ ITIH3 ELISA
 - − HCC-REAAD™ ERBB3 ELISA

REAADTM Series

- Stomach cancer detection kit
 - Early detection for Stomach cancer in blood samples
 - Incorporates the use of patented technology
 - Product brand: GC-REAAD™ ITIH3 ELISA
- Liver cancer detection kit
 - Early detection for Liver cancer in blood samples
 - Employs the use of patented technology and proprietary algorithm
 - Product brand: HCC-REAAD™ ERBB3 ELISA
- Nose cancer detection kit
 - Early detection for Nose cancer in blood samples
 - Uses proprietary technology
 - Product brand: NPC-REAAD™ EBV EA-IgA ELISA



GC-REAADTM

CE-IVD Certified in-vitro diagnostic ELISA kit for early detection of Gastric Carcinoma (stomach cancer)

GC-REAADTM - Stomach cancer detection kit

GC-REAADTM is *in-vitro* diagnostics employed the use of patented biomarker, intended for the qualitative & semi-quantitative detection of Inter-alpha-trypsin inhibitor heavy chain H3 (ITIH3) protein in human plasma for early detection of stomach cancer.

The ITIH3 Protein


- Protease inhibitor found in extracellular matrix
- One of the heavy chain belonging to family of inter-alpha-trypsin inhibitors
- Function: Ability to covalently link to hyaluronic acid (HA)
- Important factor to stabilise extracellular matrix
- Evidence showing its antagonistic relationship with tumour invasion & metastasis – tumour suppressor role

Clinical Validation for GC-REAADTM

91.0% Sensitivity

93.6% Specificity

| Validation Results | Values |
|---------------------------|--------|
| ROC Area Under Curve | 0.97 |
| Sample Size | 998 |
| Positive Cases | 144 |
| Negative Controls | 854 |
| Positive Predictive Value | 86.23% |
| Negative Predictive Value | 97.09% |





HCC-REAADTM

CE-IVD Certified in-vitro diagnostic ELISA kit for early detection of Hepatocellular Carcinoma (liver cancer)


HCC-REAAD™ - Liver cancer detection kit

HCC-REAAD™ an *in-vitro* diagnostics employed the use of patented biomarker, intended for the qualitative & semiquantitative detection of the human epidermal growth factor 3 (ERBB3) protein in human plasma. An algorithm integrating the measurement of ERBB3 and AFP for the early detection of liver cancer.


ERBB3 Protein

- A member of ERBB family of receptor tyrosine kinase (RTK)
- Plays a role in regulating cell proliferation, and differentiation
- Deregulation of ERBB signalling is associated with several human cancers

Clinical Validation for HCC-REAADTM



| Validation Results | Values |
|---------------------------|--------|
| ROC Area Under Curve | 0.976 |
| Sample Size | 379 |
| Positive Cases | 79 |
| Negative Controls | 300 |
| Positive Predictive Value | 86.81% |
| Negative Predictive Value | 95.24% |





NPC-REAADTM


CE-IVD Certified in-vitro diagnostic ELISA kit for early detection of Nasopharyngeal Carcinoma (nose cancer)

NPC-REAADTM - Nose cancer detection kit


NPC-REAAD ™ is an *in-vitro* diagnostic intended to screen for IgA antibodies specific to EBV EA as an indicator of risk for NPC

- A cocktail of 4 (out of more than 20) specific proprietary proteins derived from the EBV Early Antigen (EA) is utilised in the test kit
 - Ribonucleotide Reductase Large subunit
 - Ribonucleotide Reductase Small subunit
 - DNA Polymerase Accessory subunit
 - Novel proprietary protein
- Test results help distinguish NPC from non-NPC cases to facilitate patient diagnosis & treatment
- Clinically validated using NPC patients (instead of EBV positive patients) unlike other EBV assays

Clinical Validation for NPC-REAAD™



| Validation Results | Values |
|---------------------------|---------|
| ROC Area Under Curve | 0.99 |
| Sample Size | 143 |
| Positive Cases | 25 |
| Negative Controls | 118 |
| Positive Predictive Value | 100.00% |
| Negative Predictive Value | 98.33% |





PRODUCTS – COVID-19 TESTS



RAPIT COVID-19 IgM/IgG Test Kit


Rapid Test for Detection and Differentiation of SARS-CoV-2 IgM and IgG Antibodies

RAPIT COVID-19 IgM/IgG Test Kit


RAPIT COVID-19 IgM/IgG Test Kit is a lateral flow chromatography immunoassay intended for the qualitative detection and differentiation of IgM and IgG antibodies against SARS-CoV-2 in human plasma, serum, finger pricked whole blood or whole blood.


IgM and IgG antibodies of COVID-19 Spike and Nucleocapsid proteins are detected

RAPIT COVID-19 IgM/IgG Test Kit aids in the diagnosis of patients with suspected SARS-CoV-2 infection in conjunction with clinical presentations and the results of other laboratory tests.



Clinical Validation - RAPIT COVID-19 IgM/IgG





| Validation Results | Values |
|---------------------------|--------|
| Sample Size | 276 |
| Positive Cases | 168 |
| Negative Controls | 108 |
| Positive Predictive Value | 96.84% |
| Negative Predictive Value | 87.29% |




COVID19N-REAADTM

Enzyme-linked immunosorbent assay for qualitative detection of IgG antibodies to SARS-CoV-2 Nucleocapsid Protein

COVID19N-REAADTM


- COVID19N-REAAD™ is an enzyme—linked immunosorbent assay intended for qualitative detection of IgG antibodies to SARS-CoV-2 Nucleocapsid Protein in human serum and plasma (EDTA).
- COVID19N-REAAD™ is intended for use as an aid in identifying individuals with an adaptive immune response to SARS-CoV-2, indicating recent or prior infection.

Clinical Validation for COVID19N-REAAD™





| Validation Results | Values |
|---------------------------|---------|
| ROC Area Under Curve | 1.000 |
| Sample Size | 387 |
| Positive Cases | 56 |
| Negative Controls | 331 |
| Positive Predictive Value | 100.00% |
| Negative Predictive Value | 99.70% |





SERVICES & COLLABORATION

SERVICES

Contract Testing Services

- Offers a variety of testing services with our range of analytical instruments and expertise in immunological and molecular based assays.
- Services can be customized to support the Research & Development pipeline from Discovery to Pre-clinical phases.

Contract Development Services Contract Manufacturing Services





COLLABORATION

Licensing of IPs

Req'd for product development & commercialization

Joint product development

Evaluation studies

- Req'd for product validation
- Publication and/or report for product promotion

Supply of biological materials

- Req'd by Company for product development & manufacturing
- Req'd by Partner for R&D





Note: Services & Collaboration rendered are subject to agreement by parties involved after review over terms & mutual benefits.

BIOLOGICAL MATERIALS

Whole/Processed Blood, Plasma/Serum

- Biological and Clinical material for use in R&D/QC/Validation
- Individual clinical data
- Detailed COA & De-linked IRB Protocol
- **Anonymized Clinical Data**
- Customized requirements



| Infectious Disease/ Virology | EBV (VCA, EA, EBNA), H. Pylori, HBV DNA, Influenza A/B, Saccharomyces, Tetanus, Zika, etc. | |
|---------------------------------|---|--|
| | Hemolytic anemia, Hemophilia, Hereditary spherocytosis, Hodgkin's lymphoma, etc. | |
| Oncology | Carcinoma, Sarcoma, Lymphoma, Leukemia, etc. | |
| Others | Ulcerative Colitis, Systemic Lupus Erythematosus, etc. | |

BIOLOGICAL MATERIALS

FFPE Tissue Organs

- Database containing more than 10,000 FFPE pathological and normal tissue specimens from all organs (except brain) and diseases/cancers
- A wide range of liquid-based cytology specimens
- Raw biological pathology material or characterized by pathologists
- Sample originated from EN ISO 15189 accredited facility
- Individual clinical data
- Detailed COA & De-linked IRB Protocol
- **Anonymized Clinical Data**
- **Customized requirements**

Isolated Cell Samples from Diseased & Health Cohorts

- Immune cells (PBMCs, BMMCs, MNCs, CD4, CD25, etc.)
- Primary Respiratory cells (HuBECs, parenchymal fibroblasts, etc.)
- Primary and Immortalized cancer cells
- Joint Tissue cells (Chondrocytes, synovial fibroblasts etc.)

| Skin | Ovary | Endometrium |
|----------|------------|-------------|
| Breast | Cervical | ENT |
| Colon | Stomach | Thyroid |
| Prostate | Testicle | Lung |
| Bladder | Tissue | Liver |
| Melanoma | Pancreas | Bone Marrow |
| Kidney | Myometrium | Oesophagus |



SUMMARY

- Restalyst is an IVD company which successfully develops and produces cancer tests and infectious disease tests (COVID-19)
- Established in 2007, Restalyst is ISO 13485 certified as early as 2009
- Restalyst's IVD products are CE marked & locally registered. Products are sold globally through distributor partnerships.
- Based on the many years of experience, Restalyst offers contract services in testing, product development, production, regulatory & commercialization.
- Restalyst is also interested to collaborate with partners to supply biological materials and bring inventions from lab to be developed, produced and distributed as IVD products for the global market.

END